




nMR is **ART**



GB1 (beta1 immunoglobulin binding domain protein)

Kevin Donovan, 2016

Chemistry Department and Francis Bitter Magnet Lab
Massachusetts Institute of Technology (MIT)

Field strength: 800 MHz

Why is this your favorite spectrum?

This spectrum illustrates powerful spectroscopic techniques applied to a model sample. It was processed and contoured in a way so as not to suppress the imperfections, and is displayed with various visual filters. Presented in this manner, the spectrum is not only scientifically beautiful but rather artistically beautiful, as the imperfections add a level of abstraction by forcing us to visualize an alternate significance.